



化学薬品

# Thermo Scientific ピュラトロニック化合物 高純度無機化合物

thermo  
scientific | Chemicals

Thermo Scientific(旧AlfaAesar)輸入元・日本代理店 ヒドラス化学株式会社  
TEL.03-3258-5031 FAX.03-3258-6535 info@hydrus.co.jp www.hydrus.co.jp

# Thermo Scientific ピュラトロニック化合物

## High purity inorganics

Known worldwide as a leading manufacturer of high purity inorganic compounds, Thermo Scientific produces thousands of distinct materials to exacting standards for research, development and production applications. Custom production and packaging services are part of our regular offering. Our brands are recognized for purity and quality and are backed up by technical and sales teams dedicated to providing the best service. Our portfolio includes a selection of our wide range of high purity inorganic materials.

The Puratronic line of base metal salts, solutions and pure elements is the leading choice of pharmaceutical and high technology companies needing the basic building blocks for many research and development processes. Each Puratronic compound has a minimum purity of 99.99% (metals basis) and each is backed by a lot-specific certificate of analysis. Our Puratronic production labs are in continuous operation and can produce lots from a few grams up to hundreds of kilograms.

## Applications for inorganics

High purity products for infrared fiber optics materials for fiber optics include high-purity fluorides and sulfides. Typical purities are 99.99% to 99.995%.


Products include:

- Ammonium nitrate
- Germanium(IV) oxide
- Antimony(III) oxide
- Lead(II) fluoride
- Boric acid
- Lithium niobium oxide
- Calcium carbonate
- Lithium tetraborate
- Cesium iodide
- Tantalum(V) chloride
- Gallium(III) oxide
- Thallium(III) iodide



# Periodic table product selector:

元素をクリックして製品を参照してください。



## アルミニウム

カタログ番号	説明	パックサイズ
010622-06	Aluminum chloride hexahydrate, Puratronic, 99.9995% (metals basis)	5 g, 50 g, 250 g, 1 kg
038366-GA	Aluminum foil, 0.038 +/- 0.007mm (0.0015 +/- 0.0003in) thick, Puratronic, 99.9995% (metals basis)	50 × 200 mm, 50 × 1000 mm
011376-GH	Aluminum foil, 0.1mm (0.004in)thick, Puratronic, 99.997% (metals basis)	100 × 100 mm, 100 × 400 mm, *5 × 100 × 400 mm
043424-FI	Aluminum foil, 0.25mm (0.01in)thick, annealed, Puratronic, 99.9995% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
044331-FI	Aluminum foil, 0.25mm (0.01in)thick, annealed, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
011375-FI	Aluminum foil, 0.25mm (0.01in)thick, Puratronic, 99.997% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
043777-GH	Aluminum foil, 0.2mm (0.008in) thick, annealed, Puratronic, 99.9995% (metals basis)	100 × 100 mm, 100 × 500 mm
044330-GH	Aluminum foil, 0.2mm (0.008in)thick, annealed, Puratronic, 99.9999% (metals basis)	100 × 100 mm, 100 × 500 mm
043776-HP	Aluminum foil, 0.2mm (0.008in)thick, Puratronic, 99.997% (metals basis)	100 × 500 mm
044332-FI	Aluminum foil, 0.5mm (0.02in) thick, annealed, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
043425-GH	Aluminum foil, 0.5mm (0.02in)thick, annealed, Puratronic, 99.9995% (metals basis)	100 × 100 mm, 100 × 500 mm
011374-FI	Aluminum foil, 0.5mm (0.02in)thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
011373-FI	Aluminum foil, 1.0mm (0.04in)thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
044492-FI	Aluminum foil, 1.5mm (0.06in)thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
011372-FI	Aluminum foil, 2.0mm (0.08in) thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
040388-18	Aluminum ingot, Puratronic, 99.9995% (metals basis)	50 g, 250 g
010626-06	Aluminum nitrate hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g, 500 g
043145-14	Aluminum nitrate nonahydrate, low mercury, Puratronic, 99.999% (metals basis)	25 g, 100 g
010627-14	Aluminum oxide, Puratronic, Puratronic, 99.995% (metals basis)	25 g, 100 g, 500 g
010628-06	Aluminum phosphate, Puratronic, 99.99% (metals basis)	5 g, 25 g
010778-BQ	Aluminum rod, 22mm (0.87in) dia, Puratronic, 99.999% (metals basis)	5 cm, 25 cm
043429-BQ	Aluminum rod, 22mm (0.87in) dia, Puratronic, 99.9995% (metals basis)	5 cm, 25 cm
011493-BQ	Aluminum rod, 6.35mm (0.25in)dia, Puratronic, 99.999% (metals basis)	5 cm, *4 × cm
011478-BT	Aluminum rod, 6mm (0.2in) dia, Puratronic, 99.9965% (metals basis)	15 cm, *5 × 15 cm
010779-BQ	Aluminum rod, 9.5mm (0.37in) dia, Puratronic, 99.999% (metals basis)	5 cm, 25 cm, *4 × 25 cm
010573-18	Aluminum shot, approx. 4-8mm (0.2-0.3in), Puratronic, 99.999% (metals basis)	50 g, 1 kg
038492-14	Aluminum Silicon slug, 6.35mm (0.25in) dia × 12.7mm (0.5in) length, Puratronic, 99.999% (metals basis)	25 g
042257-09	Aluminum slug, 3.175mm (0.125in) dia × 3.175mm (0.125in) length, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
043390-09	Aluminum slug, 3.175mm (0.125in) dia × 3.175mm (0.125in) length, Puratronic, 99.9995% (metals basis)	10 g, 50 g
044326-09	Aluminum slug, 3.175mm (0.125in) dia × 3.175mm (0.125in) length, Puratronic, 99.9999% (metals basis)	10 g, 50 g
044327-09	Aluminum slug, 3.175mm (0.125in) dia × 6.35mm (0.25 in) length, Puratronic, 99.9995% (metals basis)	10 g, 50 g
042259-30	Aluminum slug, 3.175mm (0.125in) dia × 6.35mm (0.25in) length, Puratronic, 99.999% (metals basis)	250 g
043391-09	Aluminum slug, 3.175mm (0.125in) dia × 6.35mm (0.25in) length, Puratronic, 99.9995% (metals basis)	10 g
043392-09	Aluminum slug, 6.35mm (0.25in) dia × 12.7mm (0.50in) length, Puratronic, 99.9995% (metals basis)	10 g, 250 g
042324-09	Aluminum slug, 6.35mm (0.25in) dia × 12.7mm (0.50in) length, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g

## アルミニウム (続き)

カタログ番号	説明	パックサイズ
044329-09	Aluminum slug, 6.35mm (0.25in) dia x 12.7mm (0.5in) length, Puratronic, 99.9998% (metals basis)	10 g, 50 g
042323-18	Aluminum slug, 6.35mm (0.25in) dia x 6.35mm (0.25in) length, Puratronic, 99.999% (metals basis)	50 g
043357-06	Aluminum slug, 6.35mm (0.25in) dia x 6.35mm (0.25in) length, Puratronic, 99.9995% (metals basis)	5 g, 25 g
044328-14	Aluminum slug, 6.35mm (0.25in) dia x 6.35mm (0.25in) length, Puratronic, 99.9999% (metals basis)	25 g
044577-G5	Aluminum wire, 0.076mm (0.003in) dia, hard, Puratronic, 99.999% (metals basis)	5 m, 25 m
000005-G5	Aluminum wire, 0.127mm (0.005in) dia, hard, Puratronic, 99.999% (metals basis)	5 m, 25 m
044579-G5	Aluminum wire, 0.152mm (0.006in) dia, hard, Puratronic, 99.999% (metals basis)	5 m, 25 m
039730-G5	Aluminum wire, 0.18mm (0.007in) dia, Puratronic, 99.999% (metals basis)	5 m, 25 m
044578-G5	Aluminum wire, 0.1mm (0.004in) dia, hard, Puratronic, 99.999% (metals basis)	5 m, 25 m
014445-G9	Aluminum wire, 0.25mm (0.01in)dia, annealed, Puratronic, 99.999% (metals basis)	25 m, 100 m
010748-G6	Aluminum wire, 0.5mm (0.02in) dia, annealed, Puratronic, 99.999% (metals basis)	10 m, 50 m, 250 m
043426-G6	Aluminum wire, 0.5mm (0.02in)dia, annealed, Puratronic, 99.9995% (metals basis)	10 m, 50 m, 250 m
044333-G6	Aluminum wire, 0.5mm (0.02in)dia, annealed, Puratronic, 99.9999% (metals basis)	10 m
042238-G6	Aluminum wire, 0.5mm (0.02in)dia, hard, Puratronic, 99.999% (metals basis)	10 m, 50 m, 250 m
043336-H7	Aluminum wire, 0.81mm (0.032in) dia, annealed, Puratronic, 99.999% (metals basis)	250 m
043427-G5	Aluminum wire, 1.0mm (0.04in) dia, annealed, Puratronic, 99.9995% (metals basis)	5 m, 100 m
044334-G9	Aluminum wire, 1.0mm (0.04in) dia, annealed, Puratronic, 99.9999% (metals basis)	25 m
010747-G5	Aluminum wire, 1.0mm (0.04in)dia, annealed, Puratronic, 99.999% (metals basis)	5 m, 25 m, 100 m, 300 m
010765-G1	Aluminum wire, 1.5mm (0.059in)dia, annealed, Puratronic, 99.999% (metals basis)	1 m, 5 m
014441-G1	Aluminum wire, 2.0mm (0.08in)dia, annealed, Puratronic, 99.999% (metals basis)	1 m, 5 m, 25 m
043428-G1	Aluminum wire, 2.0mm (0.08in)dia, annealed, Puratronic, 99.9995% (metals basis)	1 m, 5 m, 25 m
044335-G1	Aluminum wire, 2.0mm (0.08in)dia, annealed, Puratronic, 99.9998% (metals basis)	1 m, 5 m

## アンモニウム塩

010895-09	Ammonium bis(oxalato)oxotitanate(IV) hydrate, Puratronic, 99.998% (metals basis)	1 g, 10 g, 50 g
012981-06	Ammonium bromide, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
012980-09	Ammonium carbonate, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
010632-18	Ammonium chloride, Puratronic, 99.999% (metals basis)	50 g, 250 g, 1 kg
010633-14	Ammonium dihydrogen phosphate, Puratronic, 99.995% (metals basis)	25 g, 100 g
040216-14	Ammonium fluoride hydrate, Puratronic, 99.995% (metals basis)	25 g, 100 g
043184-09	Ammonium hexafluorosilicate, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
040751-03	Ammonium iron(II) sulfate hydrate, Puratronic, 99.999% (metals basis)	1 g, 5 g, 25 g
010811-06	Ammonium molybdate (para) hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g
010846-04	Ammonium perrhenate(VII), Puratronic, 99.999% (metals basis)	2 g
010639-09	Ammonium sulfate, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
010902-06	Ammonium vanadium oxide, Puratronic, 99.995% (metals basis)	5 g, 25 g

## アンチモン

010598-09	Antimony lump, polycrystalline, Puratronic, 99.9999% (metals basis)	10 g, 50 g, 5 × 50 g, 1 kg
000428-06	Antimony shot, 1.5-5 mm (0.06-0.2 in), Puratronic, 99.9999% (metals basis)	5 g, 25 g, 100 g
011071-09	Antimony shot, 1-3 mm (0.04-0.1 in), Puratronic, 99.9999% (metals basis)	10 g, 50 g, 250 g
040372-04	Antimony(III) chloride, Puratronic, 99.999% (metals basis)	2 g
010641-18	Antimony(III) oxide, Puratronic, 99.999% (metals basis)	50 g, 100 g
010643-06	Antimony(V) oxide, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100 g
013044-03	Arsenic polycrystalline lump, 2-8 mm (0.08-0.3i n), Puratronic, 99.99999+% (metals basis)	1 g, 100 g, 1 kg
011173-04	Arsenic(III) oxide, Puratronic, 99.996% (metals basis excluding Sb)	2 g, 10 g, 50 g

## バリウム

カタログ番号	説明	パックサイズ
011126-06	Barium acetate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
010645-14	Barium carbonate, Puratronic, 99.997% (metals basis)	25 g, 100 g, 500 g
010646-14	Barium chloride dihydrate, Puratronic, 99.997% (metals basis)	25 g, 100 g, 500 g
011124-06	Barium chromate, Puratronic, 99.995% (metals basis)	5 g, 25 g
010647-14	Barium fluoride, Puratronic, 99.99% (metals basis)	25 g, 100 g, 500 g, 2.5 g
010648-06	Barium nitrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g, 500 g, 2.5 g
013009-06	Barium oxalate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
047360-03	Barium perrhenate, Puratronic, 99.99%	1 g, 5 g, 25 g
010651-06	Barium sulfate, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100 g

## ビスマス

013029-14	Bismuth germanium oxide, Puratronic, 99.9995% (metals basis)	25 g, 100 g
014442-22	Bismuth polycrystalline lump, Puratronic, 99.998% (metals basis)	100 g, 500 g, 2 kg
047229-KF	Bismuth rod, 10mm (0.39in) dia × 145mm (5.71in) long, Puratronic, 99.9999% (metals basis)	5 pc
010656-04	Bismuth(III) iodide, Puratronic, 99.999% (metals basis)	2 g, 10 g, 50 g
010657-22	Bismuth(III) nitrate hydrate, Puratronic, 99.999% (metals basis)	100 g, 500 g
011122-09	Bismuth(III) nitrate oxide, Puratronic, 99.998% (metals basis)	10 g, 50 g, 250 g
042863-18	Bismuth(III) oxide, Puratronic, 99.999% (metals basis)	50 g, 250g, 1 kg
010658-18	Bismuth(III) oxide, Puratronic, 99.9995% (metals basis)	50 g, 250 g, 1 kg

## ホウ素

010659-14	Boric acid, Puratronic, 99.9995% (metals basis)	25 g, 100 g, 500 g
011159-QE	Boron oxide, Puratronic, ( $H_2O$ 200 ppm), 99.999% (metals basis)	1 cyl, 10 cyl
011160-QE	Boron oxide, Puratronic, ( $H_2O$ 200 ppm), 99.999% (metals basis)	1 cyl, 10 cyl

## カドミウム

010660-09	Cadmium carbonate, Puratronic, 99.998% (metals basis)	10 g, 50 g
010661-09	Cadmium chloride hydrate, Puratronic, 99.998% (metals basis)	1 g, 50 g
040203-06	Cadmium fluoride, Puratronic, 99.99% (metals basis)	5 g, 25 g, 100 g
000064-FQ	Cadmium foil, 0.025 mm (0.00 1in) thick, Puratronic, 99.999% (metals basis)	25 × 100 mm, 25 × 400 mm
044127-FI	Cadmium foil, 0.05 mm (0.002 in) thick, Puratronic, 99.999% (metals basis)	50 × 50 mm
011371-FI	Cadmium foil, 0.1 mm (0.004 in) thick, Puratronic, 99.9975% (metals basis)	50 × 50 mm, 100 × 100 mm
011370-FI	Cadmium foil, 0.25 mm (0.01 in) thick, Puratronic, 99.9975% (metals basis)	50 × 50 mm, 100 × 100 mm
011369-FI	Cadmium foil, 0.5 mm (0.02 in) thick, Puratronic, 99.9975% (metals basis)	50 × 50 mm, 100 × 100 mm
011368-FI	Cadmium foil, 1.0 mm (0.04 in) thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm
011367-FF	Cadmium foil, 2.0 mm (0.08 in) thick, Puratronic, 99.998% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
010662-04	Cadmium hydroxide, Puratronic, 99.98% (metals basis)	2 g, 10 g
043442-14	Cadmium iodide, Puratronic, 99.999% (metals basis)	25 g, 100 g
010663-06	Cadmium nitrate tetrahydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
010664-14	Cadmium oxide, Puratronic, 99.998% (metals basis)	25 g, 100 g, 500 g
010120-KG	Cadmium rod, 12.7 mm (0.5 in) dia × 23 cm (9.0 in), Puratronic, 99.999+% (metals basis)	1 pc, 5 pc
043116-DM	Cadmium rod, 8 mm (0.3 in) dia, Puratronic, 99.999% (metals basis)	100 mm, 5 × 100 mm
011351-30	Cadmium shot, 1-3 mm (0.04-0.1 in), Puratronic, 99.9997% (metals basis)	250 g
010675-09	Cadmium sulfate hydrate, Puratronic, 99.996% (metals basis)	10 g, 50 g
011910-G1	Cadmium wire, 0.5 mm (0.02 in) dia, Puratronic, 99.998% (metals basis)	1 m, 5 m
000060-G4	Cadmium wire, 0.64 mm (0.025 in) dia, Puratronic, 99.999% (metals basis)	3 m, 15 m
011909-G1	Cadmium wire, 1.0 mm (0.04 in) dia, Puratronic, 99.998% (metals basis)	1 m, 5 m
011908-BU	Cadmium wire, 2.0 mm (0.08 in) dia, Puratronic, 99.998% (metals basis)	25 cm

## カルシウム

カタログ番号	説明	パックサイズ
012974-06	Calcium acetate hydrate, Puratronic, 99.9965% (metals basis)	5 g, 25 g
010679-06	Calcium carbonate, Puratronic, 99.99% (metals basis)	5 g, 25 g, 100 g, 500 g
043073-06	Calcium carbonate, Puratronic, 99.997% (metals basis)	5 g, 25 g, 100 g, 500 g
047235-06	Calcium carbonate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100g
010680-04	Calcium chloride hydrate, Puratronic, 99.9965% (metals basis)	2 g, 10 g, 50g, 250g
010681-04	Calcium fluoride, Puratronic, 99.985% (metals basis)	2 g, 10 g, 50g
011117-06	Calcium iodide hydrate, Puratronic, 99.9965% (metals basis)	5 g, 25 g
010683-06	Calcium nitrate hydrate, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100g
044515-06	Calcium nitrate hydrate, Puratronic, 99.9995% (metals basis)	5 g, 25 g, 100g
013007-06	Calcium oxalate monohydrate, Puratronic, 99.9985% (metals basis), (excluding 25 ppm max alkali earths)	5 g, 25 g
010684-04	Calcium oxide, Puratronic, 99.998% (metals basis, excluding other alkaline earth and alkali metals)	1g, 2 g, 10 g, 50g, 250g
011116-04	Calcium sulfate, anhydrous, Puratronic, 99.993% (metals basis)	2 g, 10 g, 50g

## セシウム

012117-06	Cesium carbonate, Puratronic, 99.994% (metals basis)	5 g, 25 g, 100g
012121-03	Cesium fluoride, Puratronic, 99.99% (metals basis)	1g, 10 g, 50g
011112-09	Cesium sulfate, Puratronic, 99.997% (metals basis)	10 g, 50g

## クロミウム

010151-14	Chromium pieces, electrolytic, 6 mm & down, Puratronic, 99.997% (metals basis)	25 g
010686-06	Chromium(III) chloride hexahydrate, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100g
010688-06	Chromium(III) oxide, Puratronic, 99.97% (metals basis)	5 g, 500g

## コバルト

011455-FF	Cobalt foil, 0.1mm (0.004in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011452-FI	Cobalt foil, 0.25mm (0.01in) thick, Puratronic, 99.995% (metals basis)	50 × 50 mm, 100 × 100 mm
011451-FG	Cobalt foil, 0.5mm (0.02in) thick, Puratronic, 99.995% (metals basis)	25 × 50 mm, 50 × 100 mm
011449-FF	Cobalt foil, 2.0mm (0.08in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm
010610-06	Cobalt powder, -22 mesh, Puratronic, 99.998% (metals basis)	5 g, 25 g
010867-HF	Cobalt rod, 12mm (0.47in) dia, Puratronic, 99.995% (metals basis)	50 mm, 10 mm
011445-HF	Cobalt rod, 5mm (0.2in) dia, Puratronic, 99.995% (metals basis)	50 mm, 150 mm
043227-QI	Cobalt slug, 6.35mm (0.25in) dia × 12.7mm (0.50in) length, Puratronic, 99.995% (metals basis)	1 slug, 5 slug, 25 slug
043226-QJ	Cobalt slug, 6.35mm (0.25in) dia × 6.35mm (0.25in) length, Puratronic, 99.995% (metals basis)	5 slug, 25 slug
010949-G1	Cobalt wire, 0.1mm (0.004in) dia, Puratronic, 99.995% (metals basis)	1 m, 5 m
010948-G1	Cobalt wire, 0.25mm (0.01in) dia, Puratronic, 99.995% (metals basis)	1 m, 5 m
010947-BW	Cobalt wire, 0.5mm (0.02in) dia, Puratronic, 99.995% (metals basis)	50 m
010946-BU	Cobalt wire, 1.0mm (0.04in) dia, annealed, Puratronic, 99.995% (metals basis)	25 cm, 100 cm
044525-BP	Cobalt wire, 1.5mm (0.06in) dia, annealed, Puratronic, 99.995% (metals basis)	2 cm, 10 cm, 50 m
010945-BP	Cobalt wire, 2.0mm (0.08in) dia, Puratronic, 99.995% (metals basis)	2 cm, 10 cm, 50 m
010692-09	Cobalt(II) chloride hexahydrate, Puratronic, 99.998% (metals basis)	10 g, 50 g
010693-06	Cobalt(II) fluoride tetrahydrate, Puratronic, 99.99% (metals basis)	5 g, 25 g, 50 g
010694-09	Cobalt(II) nitrate hexahydrate, Puratronic, 99.999% (metals basis)	10 g, 50 g
013012-09	Cobalt(II) oxalate dihydrate, Puratronic, 99.995% (metals basis)	10 g, 50 g, 500 g
010696-06	Cobalt(II) sulfate heptahydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g
010695-09	Cobalt(II,III) oxide, Puratronic, 99.9985% (metals basis)	10 g, 50 g

## 銅

カタログ番号	説明	パックサイズ
010950-HV	Copper foil, 0.025 mm (0.001 in) thick, Puratronic, 99.999% (metals basis)	25 × 1000 mm, 25 × 200 mm
042972-FI	Copper foil, 0.05 mm (0.002 in) thick, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm
042189-FI	Copper foil, 0.1 mm (0.004 in) thick, Puratronic, 99.999% (metals basis)	50 × 50 mm, 100 × 300 mm
042973-FI	Copper foil, 0.1 mm (0.004 in) thick, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm
000097-FY	Copper foil, 0.25 mm (0.01 in) thick, hard, Puratronic, 99.999+% (metals basis)	50 × 100 mm
011393-FI	Copper foil, 0.25 mm (0.01 in) thick, Puratronic, 99.9985% (metals basis)	50 × 50 mm, 100 × 100 mm
042974-FI	Copper foil, 0.25 mm (0.01 in) thick, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm
036395-FI	Copper foil, 0.5 mm (0.02 in) thick, Puratronic, 99.998% (metals basis), Oxygen free	50 × 50 mm, 50 × 250 mm
011392-FI	Copper foil, 0.5 mm (0.02 in) thick, Puratronic, 99.9985% (metals basis)	50 × 50 mm, 100 × 100 mm
042975-FI	Copper foil, 0.5 mm (0.02 in) thick, Puratronic, 99.9999% (metals basis)	50 × 50 mm, 100 × 100 mm
011378-FG	Copper foil, 1.0 mm (0.04 in) thick, Puratronic, 99.999% (metals basis)	25 × 50 mm, 50 × 100 mm
042976-FF	Copper foil, 1.0 mm (0.04 in) thick, Puratronic, 99.9999% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011377-FF	Copper foil, 2.0 mm (0.08 in) thick, Puratronic, 99.999% (metals basis)	25 × 25 mm, 50 × 50 mm
043822-KS	Copper plate, irregular, nominal 10 × 34 × 0.6 cm (4 × 13.5 × 0.25 in), Puratronic, 99.999% (metals basis)	1 each
012754-DJ	Copper rod, 19 mm (0.8 in) dia, Puratronic, 99.999% (metals basis)	50 mm, 250 mm, 10 mm
044695-DJ	Copper rod, 25.4 mm (1.0 in) dia, Puratronic, 99.995% (metals basis)	50 mm, 250 mm, 10 mm
000761-14	Copper rod, 3.18 mm (0.125 in) dia, random lengths, Puratronic, 99.999% (metals basis)	25 g, 100 g
011444-DO	Copper rod, 5 mm (0.2 in) dia, Puratronic, 99.999% (metals basis)	150 mm
011443-DO	Copper rod, 7 mm (0.3 in) dia, Puratronic, 99.999% (metals basis)	150 mm, *5 × 150 mm
010156-14	Copper rod, 9.5 mm (0.4 in) dia, Puratronic, 99.999% (metals basis)	25 g, 100 g, 500 g
010953-18	Copper shot, 4-6 mm (0.16-0.24 in), Puratronic, 99.999% (metals basis)	50 g, 1k g
010955-14	Copper shot, 6 mm (0.24 in) & down, Puratronic, 99.9999% (metals basis)	25 g, 100 g, 500 g
042958-09	Copper slug, 3.175 mm (0.125 in) dia × 3.175 mm (0.125 in) length, Puratronic, 99.9999% (metals basis)	10 g, 50 g
042959-09	Copper slug, 3.175 mm (0.125 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.9999% (metals basis)	10 g, 50 g
042961-22	Copper slug, 6.35 mm (0.25 in) dia × 12.7 mm (0.50 in) length, Puratronic, 99.9999% (metals basis)	100 g
042960-14	Copper slug, 6.35 mm (0.25 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.9999% (metals basis)	25 g, 100 g
010970-G6	Copper wire, 0.025 mm (0.001 in) dia, Puratronic, 99.995% (metals basis)	10 m, 50 m, 200 m
010971-G6	Copper wire, 0.1 mm (0.004 in) dia, hard, Puratronic, 99.9985% (metals basis)	10 m, 50 m
042977-G6	Copper wire, 0.1 mm (0.004 in) dia, Puratronic, 99.9999% (metals basis)	10 m, 50 m
010972-G6	Copper wire, 0.25 mm (0.01 in) dia, Puratronic, 99.9985% (metals basis)	10 m, 50 m
042978-G6	Copper wire, 0.25 mm (0.01 in) dia, Puratronic, 99.9999% (metals basis)	10 m, 50 m
010973-G5	Copper wire, 0.5 mm (0.02 in) dia, Puratronic, 99.999% (metals basis)	5 m, 25 m, 100 m
036397-BW	Copper wire, 0.5 mm (0.02 in) dia, Puratronic, 99.999% (metals basis), Oxygen Free	50 m, 250 cm
042979-G1	Copper wire, 0.5 mm (0.02 in) dia, Puratronic, 99.9999% (metals basis)	1 m, 5 m, 25 m
000098-G2	Copper wire, 0.64 mm (0.025 in) dia, Puratronic, 99.999% (metals basis), Oxygen free	2 m, 10 m
010974-G1	Copper wire, 1.0 mm (0.04 in) dia, Puratronic, 99.999% (metals basis)	1 m, 5 m, 25 m
036398-BW	Copper wire, 1.0 mm (0.04 in) dia, Puratronic, 99.999% (metals basis), Oxygen Free	50 cm, 250 cm
042980-G1	Copper wire, 1.0 mm (0.04 in) dia, Puratronic, 99.9999% (metals basis)	1 m, 5 m
044550-G1	Copper wire, 1.5 mm (0.06 in) dia, Puratronic, 99.999% (metals basis)	1 m, 5 m, 25 m
044551-G1	Copper wire, 1.5 mm (0.06 in) dia, Puratronic, 99.9999% (metals basis)	1 m, 5 m
010975-BW	Copper wire, 2.0 mm (0.08 in) dia, Puratronic, 99.999% (metals basis)	30 cm, 50 m, 150 cm, 500 cm
040752-03	Copper(I) bromide, Puratronic, 99.998% (metals basis)	1 g, 5 g
043153-06	Copper(I) iodide, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100 g
010698-06	Copper(II) chloride hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
11372098	Copper(II) fluoride hydrate, Puratronic, 99.99% (metals basis)	2 g, 10 g, 50 g

**銅 (続き)**

カタログ番号	説明	パックサイズ
010699-06	Copper(II) nitrate hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g, 500 g
010700-14	Copper(II) oxide, Puratronic, 99.995% (metals basis)	25 g, 100 g, 500 g
010701-14	Copper(II) sulfate hydrate, Puratronic, 99.999% (metals basis)	25 g, 100 g

**ガリウム**

011150-06	Gallium(III) nitrate hydrate, Puratronic, 99.999% (metals basis)	2.5 g, 5 g, 100g, 500g
042528-06	Gallium(III) oxide, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100g
011151-06	Gallium(III) oxide, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100g, 500g
011152-09	Gallium(III) sulfate hydrate, Puratronic, 99.999% (metals basis)	10 g, 50g
011144-06	Germanium pieces, 2 cm (0.8 in) & down, Puratronic, 99.9999+% (metals basis)	5 g, 25 g, 100g
011155-06	Germanium(IV) oxide, Puratronic, 99.999% (metals basis)	5 g, 25 g, 1kg

**インジウム**

011386-FI	Indium foil, 0.1 mm (0.004 in) thick, Puratronic, 99.9975% (metals basis)	50 × 50 mm, 100 × 100 mm
011385-FI	Indium foil, 0.25 mm (0.01 in) thick, Puratronic, 99.9975% (metals basis)	50 × 50 mm, 100 × 100 mm
011384-FF	Indium foil, 0.5 mm (0.02 in) thick, Puratronic, 99.9975% (metals basis)	25 × 25 mm, 50 × 50 mm, 50 × 100 mm, 100 × 200 mm
011383-FF	Indium foil, 1.0 mm (0.04 in) thick, Puratronic, 99.998% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011382-FF	Indium foil, 2.0 mm (0.08 in) thick, Puratronic, 99.998% (metals basis)	25 × 25 mm, 50 × 50 mm
010616-18	Indium ingot, polycrystalline, Puratronic, 99.9995% (metals basis)	50 g, *5 × 50 g
000151-06	Indium powder, -325 mesh, Puratronic, 99.995% (metals basis)	5 g, 25 g
040338-09	Indium shot, 1-3 mm (0.04-0.12 in), Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
011359-KM	Indium shot, 4-8 mm (0.2-0.3 in), Puratronic, 99.9999% (metals basis)	25 pc, 100 pc
010613-09	Indium shot, 5 mm (0.2 in) & down, Puratronic, 99.9995% (metals basis)	10 g, 50 g
010615-09	Indium shot, 5 mm (0.2 in) & down, Puratronic, 99.9999% (metals basis)	10 g, 50 g, 250 g
010614-09	Indium shot, 5 mm (0.2 in) & down, Puratronic, 99.9999+% (metals basis)	10 g, 50 g
040421-G3	Indium wire, 0.25 mm (0.01 in) dia, Puratronic, 99.998%+ (metals basis)	2.5 m, 10 m, 25 m
011463-G6	Indium wire, 0.5 mm (0.02 in) dia, Puratronic, 99.998%+ (metals basis)	10 m, 25 m
000155-G5	Indium wire, 0.5 mm (0.02 in) dia, Puratronic, 99.999% (metals basis)	5 m, 500 mm, 2 m
046661-G1	Indium wire, 1.0 mm (0.04 in) dia, annealed, Puratronic, 99.998%+ (metals basis)	1 m, 5 m, 25 m
011462-G1	Indium wire, 1.0 mm (0.04 in) dia, Puratronic, 99.998%+ (metals basis)	1 m, 5 m, 25 m
000470-G1	Indium wire, 1.0 mm (0.04 in) dia, Puratronic, 99.999% (metals basis)	1 m
011461-BU	Indium wire, 2.0 mm (0.08 in) dia, Puratronic, 99.998%+ (metals basis)	25 cm, 1 m, 5 m
042865-06	Indium(I) bromide, Puratronic, 99.999% (metals basis)	5 g, 25 g
010707-14	Indium(I) chloride, Puratronic, 99.999% (metals basis)	25 g
010708-06	Indium(III) nitrate hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100g, 500g
010709-09	Indium(III) oxide, Puratronic, 99.997% (metals basis)	10 g
010711-09	Indium(III) sulfate hydrate, Puratronic, 99.999% (metals basis)	10 g, 50 g, 500 g

**ヨウ素**

010619-14	Iodine, resublimed crystals, Puratronic, 99.9985% (metals basis)	25 g, 100 g, 500 g, 2.5 g
-----------	--	---------------------------

**鉄**

011381-FF	Iron foil, 0.1mm (0.004in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011380-FF	Iron foil, 0.25mm (0.01in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011379-FF	Iron foil, 0.5mm (0.02in) thick, Puratronic, 99.992% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
011448-FF	Iron foil, 1.0mm (0.04in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm, 50 × 50 mm

**鉄 (続き)**

カタログ番号	説明	パックサイズ
011447-FF	Iron foil, 2.0mm (0.08in) thick, Puratronic, 99.995% (metals basis)	25 × 25 mm, 50 × 50 mm
011447-FI	Iron foil, 2.0mm (0.08in) thick, Puratronic, 99.995% (metals basis)	50 × 50 mm
011500-DJ	Iron rod, 12mm (0.47in) dia, Puratronic, 99.995% (metals basis)	50 mm, 150 mm
011442-HT	Iron rod, 5mm (0.2in) dia, Puratronic, 99.995% (metals basis)	*5 × 150 mm
042948-G6	Iron wire, 0.075mm (0.003in) dia, Puratronic, 99.995% (metals basis)	10 m, 50 m
010938-G6	Iron wire, 0.1mm (0.004in) dia, Puratronic, 99.998% (metals basis)	10 m, 50 m
010937-G5	Iron wire, 0.25mm (0.01in) dia, Puratronic, 99.995% (metal basis)	5 m, 25 m
010936-G1	Iron wire, 0.5mm (0.02in) dia, Puratronic, 99.995% (metals basis)	1 m, 5 m
010935-BW	Iron wire, 1.0mm (0.04in) dia, Puratronic, 99.995% (metals basis)	50 cm, 250 cm
010934-BU	Iron wire, 2.0mm (0.08in) dia, Puratronic, 99.995% (metals basis)	25 cm, 100 cm
010718-14	Iron(II) oxalate dihydrate, Puratronic, 99.999% (metals basis)	25 g, 100 g, 500 g, 2.5 g
012967-06	Iron(III) sulfate heptahydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g
010715-06	Iron(III) nitrate hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g
010716-06	Iron(III) oxide, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100 g, 500 g

**鉛**

012051-FI	Lead foil, 0.1 mm (0.004 in) thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm
012050-FI	Lead foil, 0.25 mm (0.01 in) thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm
012049-FF	Lead foil, 0.5 mm (0.02 in) thick, Puratronic, 99.99% (metals basis)	25 × 25 mm, 50 × 50 mm
012048-FF	Lead foil, 1.0 mm (0.04 in) thick, Puratronic, 99.998% (metals basis)	25 × 25 mm, 50 × 50 mm
012448-FF	Lead foil, 1.0 mm (0.04 in) thick, Puratronic, 99.9995% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
012047-FF	Lead foil, 2.0 mm (0.08 in) thick, Puratronic, 99.998% (metals basis)	25 × 25 mm, 50 × 50 mm
011346-18	Lead rod, 12.7 mm (0.5 in) dia, Puratronic, 99.9998% (metals basis)	50 g, 210 g, *4 × 210 g
012450-DJ	Lead rod, 5 mm (0.2 in) dia, Puratronic, 99.9999% (metals basis)	50 mm, 250 mm
043415-15	Lead slug, 6.35 mm (0.25 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.999% (metals basis)	25 g, 100g
011902-G1	Lead wire, 0.5 mm (0.02 in) dia, Puratronic, 99.998% (metals basis)	1 m, 5 m
011901-G1	Lead wire, 1.0 mm (0.04 in) dia, Puratronic, 99.998% (metals basis)	1 m, 5 m
042861-BW	Lead wire, 2.0 mm (0.08 in) dia, Puratronic, 99.99+% (metals basis)	50 m, 250 cm
010719-06	Lead(II) acetate trihydrate, Puratronic, 99.995% (metals basis)	5 g, 50 g, 250 g, 1 kg
010720-06	Lead(II) bromide, Puratronic, 99.998% (metals basis)	5 g, 25 g, 500g
010721-06	Lead(II) carbonate, Puratronic, 99.999% (metals basis)	5 g, 25 g
010722-14	Lead(II) chloride, Puratronic, 99.999% (metals basis)	25 g, 100g
010723-14	Lead(II) fluoride, Puratronic, 99.997% (metals basis)	25 g, 100g, 2.5 g
010727-18	Lead(II) nitrate, Puratronic, 99.999% (metals basis)	50 g, 250 g, 1 kg
040753-14	Lead(II) oxide, Puratronic, 99.999% (metals basis)	25 g, 100 g, 500 g, 1 kg
010728-06	Lead(II) oxide, Puratronic, 99.9995% (metals basis)	5 g, 25 g, 100g, 500g
010730-09	Lead(II) sulfate, Puratronic, 99.999% (metals basis)	10 g, 50 g
040754-09	Lead(II) zirconium oxide, Puratronic, 99.997% (metals basis)	2 g, 10 g, 50g, 250g
010729-09	Lead(IV) oxide, Puratronic, 99.995% (metals basis)	10 g, 50 g

**リチウム**

012982-06	Lithium acetate hydrate, Puratronic, 99.998% (metals basis)	5 g, 25 g
010733-14	Lithium bromide hydrate, Puratronic, 99.995% (metals basis)	25 g, 100 g
010734-18	Lithium carbonate, Puratronic, 99.998% (metals basis excluding Ca)	50 g, 250 g, 1 kg
010735-14	Lithium chloride hydrate, Puratronic, 99.996% (metals basis)	25 g, 100 g
010736-09	Lithium fluoride, Puratronic, 99.99% (metals basis)	1g, 10 g, 50 g, 250 g
010738-06	Lithium iodide hydrate, Puratronic, 99.995% (metals basis)	5 g, 25 g
010739-14	Lithium metaborate, Puratronic, 99.997% (metals basis)	25 g, 100 g
010741-14	Lithium niobium oxide, Puratronic, 99.998% (metals basis excluding Ta), Ta <50ppm	25 g, 100 g, 500 g

## リチウム (続き)

カタログ番号	説明	パックサイズ
010742-09	Lithium nitrate hydrate, Puratronic, 99.999% (metals basis)	10 g, 50 g
010743-14	Lithium phosphate, Puratronic, 99.99% (metals basis)	25 g, 100 g
012983-06	Lithium sulfate hydrate, Puratronic, 99.996% (metals basis)	5 g, 25 g
010744-14	Lithium tetraborate, Puratronic, 99.998% (metals basis)	25 g, 1 kg

## マグネシウム

010795-14	Magnesium acetate tetrahydrate, Puratronic, 99.997% (metals basis)	25 g, 100 g
010796-06	Magnesium carbonate hydroxidehydrate, Puratronic, 99.996% (metals basis)	5 g, 100 g
010797-14	Magnesium chloride hexahydrate, Puratronic, 99.999% (metals basis)	25 g, 100 g
040755-09	Magnesium fluoride, Puratronic, 99.99% (metals basis)	10 g, 50 g
010799-06	Magnesium nitrate hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g, 500 g, 2.5 g
010800-09	Magnesium oxide, Puratronic, 99.995% (metals basis)	10 g, 50 g, 250 g
010801-09	Magnesium sulfate hydrate, Puratronic, 99.997% (metals basis)	10 g, 50 g
010630-06	Magnesium turnings, 1 cm (0.4 in) & down, Puratronic, 99.98% (metals basis)	5 g, 25 g, 100 g, 500 g

## マンガン

010802-06	Manganese(II) acetate tetrahydrate, Puratronic, 99.999% (metals basis)	5 g
040756-06	Manganese(II) carbonate, Puratronic, 99.985% (metals basis)	5 g, 25, 100g
010804-09	Manganese(III) chloride tetrahydrate, Puratronic, 99.999% (metals basis)	10 g, 50g
047109-14	Manganese(III) fluoride, Puratronic, 99.99% (metals basis)	25 g, 100g
010806-04	Manganese(III) nitrate hydrate,Puratronic, 99.995% (metals basis)	2 g, 10 g, 50g
010807-06	Manganese(III) sulfate hydrate,Puratronic, 99.999% (metals basis)	5 g, 25 g
010805-06	Manganese(IV) oxide, Puratronic, 99.996% (metals basis)	5 g, 25 g, 100g

## 水銀

010808-14	Mercury(II) chloride, Puratronic, 99.999% (metals basis)	25 g, 100g, 500g
012992-06	Mercury(II) sulfide, black, Puratronic, 99.999% (metals basis)	5 g, 25 g

## モリブデン

012972-22	Molybdenum powder, -22 mesh, Puratronic, 99.9975% (metals basis exluding W), W<1000ppm	100 g, 2.5 g
039686-14	Molybdenum powder, -325 mesh, Puratronic, 99.997% (metals basis excluding W), W<130ppm	25 g, 100 g
011475-DO	Molybdenum rod, 5mm (0.2in) dia, Puratronic, 99.97% (metals basis)	150 mm, *10 × 150 mm
010812-14	Molybdenum(VI) oxide, Puratronic, 99.998% (metals basis excluding W)	25 g, 100 g
012930-14	Molybdenum(VI) oxide, Puratronic, 99.9995% (metals basis excluding W)	25 g, 500 g

## ニッケル

012046-FI	Nickel foil, 0.1 mm (0.004 in) thick, Puratronic, 99.994% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
012045-FG	Nickel foil, 0.25 mm (0.01 in) thick, Puratronic, 99.994% (metals basis)	25 × 50 mm, 50 × 100 mm, 100 × 250 mm
012044-FG	Nickel foil, 0.5 mm (0.02 in) thick, Puratronic, 99.994% (metals basis)	25 × 25 mm, 100 × 100 mm
012043-FG	Nickel foil, 1.0 mm (0.04 in) thick, Puratronic, 99.9945% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
012042-FG	Nickel foil, 2.0 mm (0.08 in) thick, Puratronic, 99.9945% (metals basis)	25 × 25 mm, 50 × 50 mm
010674-22	Nickel powder, -120 mesh, Puratronic, 99.996% (metals basis)	100 g, 500 g, 2.5 g, 1 kg
012966-14	Nickel powder, low carbon, Puratronic, 99.999% (metals basis), C<100ppm	25 g, 100 g, 1 kg
011499-DJ	Nickel rod, 12 mm (0.47 in) dia,Puratronic, 99.995% (metals basis)	50 mm, 150 mm, 10 mm
042835-DJ	Nickel rod, 5 mm (0.2 in) dia, Puratronic, 99.999% (metals basis)	50 mm, 150 mm
042333-09	Nickel slug, 3.175 mm (0.125 in) dia × 3.175 mm (0.125 in) length, Puratronic, 99.995% (metals basis)	10 g, 50 g, 250 g
042332-18	Nickel slug, 3.175 mm (0.125 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.995% (metals basis)	50 g

## ニッケル (続き)

カタログ番号	説明	パックサイズ
042330-QI	Nickel slug, 6.35 mm (0.25 in) dia x 12.7 mm (0.50 in) length, Puratronic, 99.995% (metals basis)	1 slug, 5 slug, 25 slug
042331-QI	Nickel slug, 6.35 mm (0.25 in) dia x 6.35 mm (0.25 in) length, Puratronic, 99.995% (metals basis)	1 slug
010933-G5	Nickel wire, 0.05 mm (0.002 in) dia, Puratronic, 99.99% (metals basis)	5 m, 25 m
010932-G6	Nickel wire, 0.1 mm (0.004 in) dia, Puratronic, 99.994% (metals basis)	10 m, 50 m
010931-G2	Nickel wire, 0.25 mm (0.01 in) dia, Puratronic, 99.994% (metals basis)	2 m, 10 m, 50 m
010930-G1	Nickel wire, 0.5 mm (0.02 in) dia, Puratronic, 99.995% (metals basis)	1 m, 5 m
043761-G5	Nickel wire, 0.5 mm (0.02 in) dia, Puratronic, 99.999% (metals basis)	5 m
010929-BU	Nickel wire, 1.0 mm (0.04 in) dia, annealed, Puratronic, 99.995% (metals basis)	25 cm, 100 cm, *5 x 100 cm
043762-BU	Nickel wire, 1.0 mm (0.04 in) dia, Puratronic, 99.999% (metals basis)	25 cm, 500 cm
010928-BQ	Nickel wire, 2.0 mm (0.08 in) dia, Puratronic, 99.995% (metals basis)	5 cm, 25 cm, 100 cm
043763-BQ	Nickel wire, 2.0 mm (0.08 in) dia, Puratronic, 99.999% (metals basis)	5 cm, 25 cm, 100 cm
010813-06	Nickel(II) acetate tetrahydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g
043185-14	Nickel(II) chloride hexahydrate, Puratronic, 99.9995% (metals basis)	25 g, 100 g
010815-06	Nickel(II) chloride hydrate, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100 g
010816-06	Nickel(II) nitrate hexahydrate, Puratronic, 99.9985% (metals basis)	5 g, 250 g
040757-06	Nickel(II) oxalate hydrate, Puratronic, 99.9985% (metals basis)	5 g, 25 g
010819-06	Nickel(II) oxide, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100 g
010820-14	Nickel(II) sulfate hydrate, Puratronic, 99.9985% (metals basis)	25 g, 100 g

## ニオブル

040510-09	Niobium powder, -325 mesh, Puratronic, 99.99% (metals basis excluding Ta)	10 g, 50 g, 250 g
035365-QI	Niobium slug, 6.35 mm (0.25 in) dia x 12.7 mm (0.50 in) length, Puratronic 99.99% (metals basis excluding Ta), Ta 300ppm	1 slug
043729-06	Niobium(V) chloride, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
010822-18	Niobium(V) oxide, Puratronic, 99.9985% (metals basis)	50 g, 250 g, 1 kg

## リン

010670-03	Phosphorus lump, red, Puratronic, 99.999% (metals basis)	1g, 10 g, 50 g, 500 g
000248-03	Phosphorus lump, red, Puratronic, 99.999+% (metals basis)	1g, 10 g, 50 g

## カリウム

012975-06	Potassium acetate hydrate, Puratronic, 99.997% (metals basis)	5 g, 25 g
042898-09	Potassium bis(oxalato)oxotitanate(IV) hydrate, Puratronic, 99.997% (metals basis)	2 g, 10 g, 50 g
010837-09	Potassium bromide, Puratronic, 99.999% (metals basis)	10 g, 50 g
010838-06	Potassium carbonate, Puratronic, 99.997% (metals basis)	5 g, 25 g, 100g, 2.5 g
010839-18	Potassium chloride, Puratronic, 99.997% (metals basis)	50 g, 250 g
010840-14	Potassium dihydrogen phosphate, Puratronic, 99.999% (metals basis)	25 g, 100 g
010842-06	Potassium iodide, Puratronic, 99.998% (metals basis)	5 g, 25 g, 100 g
010844-06	Potassium nitrate, Puratronic, 99.995% (metals basis)	5 g, 50 g, 250 g
010845-06	Potassium sulfate, Puratronic, 99.997% (metals basis)	5 g, 25 g

## レニウム

010673-03	Rhenium powder, -22 mesh, Puratronic, 99.997% (metals basis)	1 g, 5 g
047316-06	Rhenium(III) chloride, Puratronic, 99.99%	5 g, 250 mg

## ルビジウム

カタログ番号	説明	パックサイズ
011105-03	Rubidium carbonate, Puratronic, 99.975% (metals basis)	1 g, 25 g
010851-09	Rubidium chloride, Puratronic, 99.975% (metals basis)	10 g, 50 g
012873-09	Rubidium fluoride, Puratronic, 99.975% (metals basis)	2 g, 10 g
010852-03	Rubidium nitrate, Puratronic, 99.975% (metals basis)	1 g, 10 g, 50 g

## セレニウム

000290-18	Selenium shot, amorphous, 1-6 mm (0.04-0.24 in), Puratronic, 99.999+% (metals basis)	50 g, 250 g
010603-18	Selenium shot, amorphous, 2-6 mm (0.08-0.2 in), Puratronic, 99.999% (metals basis)	50 g, 250 g, 1 kg
010854-06	Selenium(IV) oxide, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
010855-14	Selenous acid, Puratronic, 99.999% (metals basis)	25 g, 100 g

## ケイ素

010856-14	Silicon(IV) oxide, Puratronic, 99.999% (metals basis)	25 g, 100 g, 1 kg
-----------	---	-------------------

## ストロンチウム

012990-06	Strontium acetate hydrate, Puratronic, 99.9965% (metals basis)	5 g, 25 g
036618-14	Strontium carbonate, Puratronic, 99.994% (metals basis)	25 g
010877-06	Strontium chloride hexahydrate, Puratronic, 99.9965% (metals basis)	5 g, 25 g
010878-06	Strontium fluoride, Puratronic, 99.99% (metals basis)	0.5 g, 5 g, 25 g, 100g
010879-06	Strontium nitrate, Puratronic, 99.9965% (metals basis)	5 g, 25 g, 100g, 2.5 g
043193-14	Strontium nitrate, Puratronic, 99.999% (metals basis)	25 g, 100 g
012870-06	Strontium sulfate, Puratronic, 99.9965% (metals basis)	5 g, 25 g

## 硫黄

010755-09	Sulfur pieces, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
010343-14	Sulfur pieces, Puratronic, 99.9995% (metals basis)	25 g, 100 g, 500 g
010874-14	Sodium sulfate, Puratronic, 99.9955% (metals basis)	25 g, 100 g

## タンタル

010344-14	Tantalum powder, -100 mesh, Puratronic, 99.98% (metals basis)	25 g, 100 g
010753-06	Tantalum powder, -22 mesh, Puratronic, 99.95% (metals basis)	5 g, 25 g, 100g, 500g
010345-14	Tantalum powder, -325 mesh, Puratronic, 99.97% (metals basis), Nb 50ppm	25 g, 100 g, 500 g, 2.5 g
011441-DM	Tantalum rod, 4.5 mm (0.18 in) dia, Puratronic, 99.96% (metals basis), NB nominal	100 mm
043731-06	Tantalum(V) chloride, anhydrous, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100 g
014708-14	Tantalum(V) chloride, Puratronic, 99.99% (metals basis)	25 g, 100 g
010880-06	Tantalum(V) chloride, Puratronic, 99.99% (metals basis)	5 g, 25 g, 100 g
010881-14	Tantalum(V) oxide, Puratronic, 99.993% (metals basis excluding Nb)	25 g

## テルル

012607-14	Tellurium broken ingot, Puratronic, 99.9999+% (metals basis)	25 g, 100 g, 500 g
043960-18	Tellurium powder, -18+60 mesh, Puratronic, 99.999% (metals basis)	50 g, 250 g
010882-14	Tellurium(IV) oxide, Puratronic, 99.9995% (metals basis)	25 g, 500 g

## タリウム

043001-09	Thallium(I) iodide, Puratronic, 99.9995% (metals basis)	10 g, 50 g
-----------	---	------------

## スズ

011460-FI	Tin foil, 0.1 mm (0.004 in) thick, Puratronic, 99.998% (metals basis)	100 × 100 mm, 100 × 500mm
011459-FI	Tin foil, 0.25 mm (0.01 in) thick, Puratronic, 99.998% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 500 mm
011458-GH	Tin foil, 0.5 mm (0.02 in) thick, Puratronic, 99.9985% (metals basis)	100 × 100 mm

## スズ(続き)

カタログ番号	説明	パックサイズ
011457-FG	Tin foil, 1.0 mm (0.04 in) thick, Puratronic, 99.9985% (metals basis)	25 × 50 mm, 50 × 100 mm, 100 × 250 mm
011456-FG	Tin foil, 2.0 mm (0.08 in) thick, Puratronic, 99.9985% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
045481-06	Tin pieces, Puratronic, 99.9999+% (metals basis)	5 g, 25 g, 100g
000788-06	Tin rod, 13 mm (0.51 in) dia, Puratronic, 99.998+% (metals basis)	5 g, 25 g
011472-DP	Tin rod, 6 mm (0.2 in) dia, Puratronic, 99.9985% (metals basis)	160 mm, *5 × 160 mm
011009-14	Tin shot, 3 mm (0.1 in), Puratronic, 99.9999% (metals basis)	25 g, 100 g
014134-KK	Tin slug, 11 mm (0.43 in) dia × 4 mm (0.2 in) thick, Puratronic, 99.999% (metals basis)	1 pc, 50 pc
043416-14	Tin slug, 6.35 mm (0.25 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.995% (metals basis)	25 g, 100 g
011907-G5	Tin wire, 0.25 mm (0.01 in) dia, Puratronic, 99.998% (metals basis)	5 m, 25 m
011906-G1	Tin wire, 0.5 mm (0.02 in) dia, Puratronic, 99.9985% (metals basis)	1 m, 5 m, 25 m
011905-BW	Tin wire, 1.0 mm (0.04 in) dia, Puratronic, 99.998% (metals basis)	50 m, 250 cm, 1000 cm
011904-BU	Tin wire, 2.0 mm (0.08 in) dia, Puratronic, 99.9985% (metals basis)	25 cm, 100 cm, 500 cm
010894-06	Tin(II) chloride hydrate, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100 g
010891-06	Tin(IV) oxide, Puratronic, 99.996% (metals basis)	5 g, 25 g

## チタン

036680-KG	Titanium plate, 6mm (0.2in) thick, 100x100mm (3.9x3.9in), Puratronic, 99.999% (metals basis)	1 pc
036681-BQ	Titanium rod, 12.7mm (0.5in) dia, Puratronic, 99.999% (metals basis)	5 cm, 25 cm
042394-06	Titanium slug, 3.175mm (0.125in) dia × 3.175mm (0.125in) length, Puratronic, 99.995% (metals basis)	5 g
042393-06	Titanium slug, 3.175mm (0.125in) dia × 6.35mm (0.25in) length, Puratronic, 99.995% (metals basis)	5 g, 25 g
042391-QJ	Titanium slug, 6.35mm (0.25in) dia × 12.7mm (0.50in) length, Puratronic, 99.995% (metals basis)	5 slug, 25 slug
042392-QJ	Titanium slug, 6.35mm (0.25in) dia × 6.35mm (0.25in) length, Puratronic, 99.995% (metals basis)	5 slug, 25 slug
010897-09	Titanium(IV) oxide, Puratronic, 99.995% (metals basis)	10g, 50g, 250g

## タンクステン

039691-04	Tungsten powder, APS 4-6 micron, Puratronic, 99.999% (metals basis), Mo <1ppm	2 g, 10 g, 50 g
013398-14	Tungsten(VI) oxide, Puratronic, 99.998% (metals basis excluding Mo), Mo 100ppm	25 g, 100 g, 1 kg

## 亜鉛

012960-06	Zinc bromide hydrate, Puratronic, 99.999% (metals basis)	5 g, 25 g, 100 g
035782-09	Zinc chloride hydrate, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g
012986-06	Zinc fluoride, Puratronic, 99.995% (metals basis)	5 g, 25 g, 100 g
011912-FI	Zinc foil, 0.1 mm (0.004 in) thick, Puratronic, 99.994% (metals basis)	50 × 50 mm, 100 × 100 mm, 100 × 300 mm
011913-FI	Zinc foil, 0.25 mm (0.01 in) thick, Puratronic, 99.994% (metals basis)	50 × 50 mm, 100 × 100 mm
011914-GH	Zinc foil, 0.5 mm (0.02 in) thick, Puratronic, 99.994% (metals basis)	100 × 100 mm
011915-FF	Zinc foil, 1.0 mm (0.04 in) thick, Puratronic, 99.9985% (metals basis)	25 × 25 mm, 50 × 50 mm, 100 × 100 mm
043504-FF	Zinc foil, 2 mm (0.08 in) thick, Puratronic, 99.999% (metals basis)	25 × 25 mm, 50 × 50 mm
011136-18	Zinc nitrate hexahydrate, Puratronic, 99.998% (metals basis)	50 g, 250 g
012985-14	Zinc oxalate dihydrate, Puratronic, 99.999% (metals basis)	25 g, 100 g
011137-09	Zinc oxide, Puratronic, 99.999% (metals basis)	10 g, 50 g, 250 g, 500 g
012855-14	Zinc oxide, Puratronic, 99.9995% (metals basis)	25 g, 100 g
013013-06	Zinc phosphate (ortho), Puratronic, 99.995% (metals basis)	5 g, 25 g
010760-18	Zinc shot, 1-6 mm (0.04-0.24 in), Puratronic, 99.9999% (metals basis)	50 g, 250 g, 1 kg

**亜鉛 (続き)**

カタログ番号	説明	パックサイズ
043417-14	Zinc slug, 6.35 mm (0.25 in) dia × 6.35 mm (0.25 in) length, Puratronic, 99.993% (metals basis)	25 g, 100 g
011138-09	Zinc sulfate hydrate, Puratronic, 99.999% (metals basis)	10 g, 50 g
012055-G2	Zinc wire, 0.25 mm (0.01 in) dia, Puratronic, 99.994% (metals basis)	2 m, 10 m, 50 m
012054-G1	Zinc wire, 0.5 mm (0.02 in) dia, Puratronic, 99.994% (metals basis)	1 m, 5 m, 25 m
012053-BW	Zinc wire, 1.0 mm (0.04 in) dia, Puratronic, 99.9985% (metals basis)	50 m, 250 cm, *2 × 500 cm
042637-BU	Zinc wire, 2.0 mm (0.08 in) dia, Puratronic, 99.999% (metals basis)	25 cm, 100 cm

**ジルコニウム**

011141-22	Zirconium(IV) oxide, Puratronic, 99.978% (metals basis)	100 g, 500 g
-----------	---	--------------

**Related resources**

**Analytical standards**  
for trace elemental analysis



**High purity acids**  
for trace elemental analysis



**Materials for electronics and semiconductor industries**



**Energy storage and battery systems**

■ 詳しい情報は [thermofisher.com/chemicals](http://thermofisher.com/chemicals) をご覧いただけ、  
[chemicals.japan@thermofisher.com](mailto:chemicals.japan@thermofisher.com) までメールでお問い合わせください。

研究用にのみ使用できます。診断用には使用いただけません。実際の価格は、弊社販売代理店までお問い合わせください。価格、製品の仕様、外観、記載内容は予告なしに変更する場合がありますのであらかじめご了承ください。© 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. BN018JA0623